



STORAGE OF FLAMMABLE PRODUCTS



AC 1200 S



AC 600 S

SMART LINE



EN 14470-1
EN 16121
EN 16122



- Completely produced in pickled steel, 10/10 mm thickness, painted with epoxy acid-proof powders. Thickness of painting from a minimum of 60µm (micron) to a maximum of 250µm (micron).
- External monolithic construction, internal body resistant to corrosive products.
- Innovative fire-proof insulation made of highly fire resistant material, completely ecological.
- Double, inflatable and anti-dust trimmings, and with protection against "cold and hot fumes".
- Door's closure with thermic device at 50° C. Not electric.
- Mechanical certified safety valves with automatic closure at 70 ± 10° C.
- Tray-shaped shelves adjustable in height by racks, and bottom basin, all made of electrogalvanized steel painted with epoxy acid-proof powders.
- Key lock with extractable cylinder (on demand "passe-partout kit" master key in case of lost keys)
- Internal patented safety catch ignis-lock®.
- "Product and public liability insurance" valid for all lifelong of the cabinet.

Product Code	External dimensions LxDxH	Type	Weight	Capacity liters	No. of bottles for each shelf	Shelves	Door locking system	Palletized and certified base
AC 1200 S	1200x640x1950 mm	90	380 kg	160 lt	21x2,5 lt	3	Yes	Yes
AC 900 S	900x640x1950 mm	90	305 kg	113 lt	15x2,5 lt	3	Yes	Yes
AC 600 S	680x640x1950 mm	90	252 kg	100 lt	10x2,5 lt	3	Yes	Yes
AC 1200/130 S	1200x640x1325 mm	90	290 kg	120 lt	40x1 lt	3	Yes	Yes
AC 900/130 S	900x640x1325 mm	90	218 kg	96 lt	32x1 lt	3	Yes	Yes



AC 1200 CM



AC 600 CM

CM LINE CME LINE T30 LINE



EN 14470-1
EN 16121
EN 16122

- Completely produced in pickled steel, 10/10 mm thickness, painted with epoxy acid-proof powders. Thickness of painting from a minimum of 60µm (micron) to a maximum of 250µm (micron).
- External monolithic construction, internal body resistant to corrosive products.
- Innovative fire-proof insulation made of highly fire resistant material, completely ecological.
- Double, inflatable and anti-dust trimmings, and with protection against "cold and hot fumes".
- Door's closure with thermic device at 50° C. Not electric.
- Mechanical certified safety valves with automatic closure at 70 ± 10° C.
- Tray-shaped shelves adjustable in height by racks (depending on the model), and bottom basin, all made of electrogalvanized steel painted with epoxy acid-proof powders.
- Key lock with extractable cylinder (on demand "passe-partout kit" master key in case of lost keys).
- Patented safety catch ignis-lock®.
- "Product and public liability insurance" valid for all lifelong of the cabinet.

Product Code	External dimensions LxDxH	Type	Weight	Capacity liters	No. of bottles for each shelf	Shelves	Door locking system	Palletized and certified base
AC 1200 CM	1180x640x1950 mm	90	490 kg	160 lt	21x2,5 lt	3	Yes	Yes
AC 600 CM	680x640x1950 mm	90	350 kg	100 lt	9x2,5 lt	3	Yes	Yes
AC 1200 T30	1180x640x1950 mm	30	415 kg	160 lt	21x2,5 lt	3	Yes	Yes
AC 600 T30	680x640x1950 mm	30	350 kg	100 lt	9x2,5 lt	3	Yes	Yes
AC 1200 CM E	1180x640x1950 mm	90	490 kg	160 lt	21x2,5 lt	3	No	Yes
AC 600 CM E	680x640x1950 mm	90	350 kg	100 lt	9x2,5 lt	3	No	Yes
AC 600/50 CM	595x510x600 mm	90	100 kg	15 lt	6x2,5 lt	1	Yes	-
AC 600/50 CMD	595x510x600 mm	90	110 kg	15 lt	6x2,5 lt	1 drawer	-	-
AC 900/50 CM	895x510x600 mm	90	145 kg	25 lt	10x2,5 lt	1	Yes	-
AC 900/50 CMD	895x510x600 mm	90	146 kg	20 lt	8x2,5 lt	1 drawer	-	-
AC 1200/50 CM	1100x510x600 mm	90	170 kg	30 lt	12x2,5 lt	1	Yes	-
AC 1200/50 CMD	1100x510x600 mm	90	170 kg	30 lt	12x2,5 lt	1 drawer	-	-
AC 1200/50 CMDD	1100x510x600 mm	90	170 kg	30 lt	6x2,5 lt	2 drawers	-	-
MINIBOX	400x400x445 mm	30	35 kg	10 lt	4x2,5 lt	1	-	-

(*) Not provided on cabinets with drawer.

SAFETYBOX® FLAMMABLES

Accessories

PRODUCT	FLAMMABLES	CHEMICALS	KEMFIRE	COMPRESSED GAS
1 KIT05 - CERTIFIED DEVICE TO POUR LIQUIDS	•		•	
2 UNICUM ASPIRATOR AND FILTER	•			•
3 EXHAUST FUMES CONVEYOR	•			•
4 ANTI-SPARK OR ATEX ASPIRATORS	•		•	•
5 INTERNAL FILTRATION SYSTEM KIT F	•		•	•
6 CUBETTO® CONTROLLED COOLING UNIT	•		•	•
7 ECOBOX® PVC OR INOX CONTAINER FOR POLLUTING LIQUIDS	•	•	•	
8 ADDITIONAL SHELVES IN STAINLESS STEEL	•	•	•	
9 ADDITIONAL SHELVES IN POLYPROPYLENE	•	•	•	
10 EASY-GOING® ROLLER SYSTEM FOR SHORT MOVING	•			
11 EASY ROLLER® ROLLING SYSTEM	•			
12 ACTIVATED CHARCOAL FILTER		•	•	
13 KIT01 FAN AND FILTER FOR CABINETS			•	
14 PASSE-PARTOUT IN CASE OF LOST KEYS	•		•	
15 ANTIMICROBIAL PAINT		•	•	
16 E-FLOW® ASPIRATION SYSTEM		•	•	
17 PARTITION WALL KIT		•	•	
18 CONTAINER FOR SMALL GAS CYLINDERS				•



Via E. Majorana, 12
31025 Santa Lucia di Piave - TV - Italy
Tel. +39 0438 451054 - Fax +39 0438 455784
info@safetybox.it - www.safetybox.it



Numero Verde
800 604020

The manufacturer reserves the right to modify the product without any notice. Reproduction in whole or in part of the contents of this catalog is forbidden.

LABOR SECURITY SYSTEM



Labor Security System is leading company in the planning and production of safety cabinets SAFETYBOX® for the storage of flammable products, compressed gas cylinders, acids and bases, as well as charcoal filters CARBOX®. Pioneering, innovation and constant attention to customer's needs are the features that distinguish the company for over 25 years.



MISSION

We guarantee the highest standards of quality for customers with a reliable and CERTIFIED as well as functional product, technologically competitive, in compliance with the most stringent European standards for safety.

CERTIFICATIONS

Since 1995 the company is certified UNI EN ISO 9001 by TÜV SÜD-Italy, known internationally for the rigidity of the regulations which both company organization and product realization should meet.

SAFETYBOX® products are certified by TÜV Production monitored and Safety tested, which guarantees both the quality of the product and the control of the production phase: from materials to the final testing.

Safety cabinets are CERTIFIED in compliance with EN 14470-1, EN 14470-2, EN 16121, EN 16122 and EN 61010-1 norms. The models for flammable storage and gas cylinders lines have all been tested for fire resistance in Ministry licensed research institutes.



CUSTOMER FOCUS AND RESEARCH

Innovation, research and tradition reliability are recurring themes, consolidated over 25 years of activity. The utmost attention to the customer is manifested in the study and production of products always in line with the needs of the customer, in compliance with current safety regulations, providing a pre- and post-sales assistance always up-to-date, competent, precise and prompt.

The Company provides the proper training and information to resellers through in-depth courses that take place at Labor Security System showroom. Currently present with its dealers across the country and in continuous expansion in the international market, the Company has a large warehouse in which over 400 different cabinets are stored, ready to be shipped to the customer within 24 hours from receipt of order.

