



MINERAL  
INNOVATIVE  
TECHNOLOGIES

PRODUCTIVITY IN MINERALS

# KETTLE FURNACE

**MIT-GROUP**  
MINERAL INNOVATIVE TECHNOLOGIES GROUP



# KETTLE FURNACE

Top loading benchtop smelter. Supplied standard with one crucible and lifting / pouring tong.

## OPERATING TEMPERATURE

- Max Temperature: 1150°C

## INSTRUMENTATION

- Digital temperature controller
- Type (k) thermocouple
- On & Off switch.

## DIMENSIONS

- 275mm(l) x 275mm(w) x 430mm(h)





## CONTACT US

Mineral Innovative Technologies (PTY) Ltd

**Tel:** +27 861 500 108

**Email:** [info@m-i-t.co.za](mailto:info@m-i-t.co.za)

**Web:** [www.m-i-t.co.za](http://www.m-i-t.co.za)  
[www.mit-group.co.za](http://www.mit-group.co.za)

